



**Klausz Dental
Laboratories Ltd.**

Working **harder and smarter** for your practice!

www.klauszdentallab.com • klauszdl@uniserve.com • Tel: 416-630-6262 • 1-800-410-4392
40 Kodiak Crescent • Units 11-12 • Toronto, ON • Canada M3J 3G5

An O.I.C.S. Specialty Laboratory
Occlusion - Implants - Cosmetic- Sleep

ACU-Flow Bite Registration Material Now Available from Klausz Dental Laboratory

We are always looking for ways we can improve on the quality of the restorations and services we offer our clients. When we find a material, process, or technique that can help us deliver a better result to your patients, we're eager to share what we have learned.

We are proud to announce you can now purchase ACU-Flow Bite Registration Material directly from Klausz Dental Laboratories.

The advantages of this material are clear:

- A High Durometer Product
- Will Not Slump
- Will Not Fragment
- Easy to Trim
- Will Not Distort under Pressure of Articulation
- Working time - 15 sec.
- Minimum intraoral setting time - 45 sec.
- Maximum total cure time - 1 minute

A Syringeable vinyl poly siloxane material for achieving highly accurate and detailed bite registrations.

It's the perfect companion product to use with the Aqualizer Bite Registration Technique, or any other Functional or Centric Relation Bite Registration Technique.

It's the material of choice of many of today's leading clinicians.



We've worked hard to earn your trust, and it's a trust we do not take lightly. Every day we strive to exceed your expectations and being able to provide you with another great product is just another way we are,

“Working Harder and Smarter for your Practice!”