



**Klausz Dental
Laboratories Ltd.**

Working **harder and smarter** for your practice!

www.klauszdentalab.com • klauszdl@uniserve.com • Tel: 416-630-6262 • 1-800-410-4392
40 Kodiak Crescent • Units 11-12 • Toronto, ON • Canada M3J 3G5



An O.I.C.S. Specialty Laboratory *Occlusion - Implants - Cosmetics - Sleep*

Dental Orthotics To Meet Your Specific Needs

Through constant innovation and education, Klausz Dental Laboratories has grown into an Occlusion, Implant, Cosmetic, and Sleep Specialty Laboratory, and we have solidified our position as a leader in these fields through various certifications. We are the only dental laboratory that builds advanced occlusal concepts into every restoration, orthotic, and dental appliance we custom fabricate.

One of the fastest growing segments of study in dentistry today is in the treatment of T.M.J. Disorders, Orofacial Pain and Sleep Disorders. Klausz Dental Laboratories is proud of our involvement in providing dentists with guidance and educational opportunities in this developing field. We can custom fabricate a wide variety of dental orthotics, from nighttime anterior deprogrammers to daytime orthotics, following the protocol set out by many of the leading authorities in dentistry today.

To make it easier for you to communicate and prescribe the latest in orthotics and dental appliances, we have created a new Orthotics Reference Page. Laminated for easy cleaning, feel free to request multiple copies, so they can be kept in each operatory for easy reference. Along with the Orthotics Reference Page is a new corresponding Orthotics Prescription Form. We are pleased to be able to provide additional copies, at your request.



Mandibular Day Orthotics

Mandibular Positioning Appliance

- * Custom PMT Base or Flexible Thermoplastic Base
- * Optional ball clasp 3-4 per side
- * Full coverage over Mid posterior teeth
- * Optional anterior wrap-over incisal coverage
- * Lingual anatomical ridge 30% of Mx post cusp tips index
- * No contact of Mx buccal cusps
- * No protrusive and non-working interferences

Maxillary Night Orthotics

Anterior Deprogrammer

- * Custom PMT Base or Flexible Thermoplastic Base
- * Optional with ball clasp between 5-6's
- * Anterior acrylic pad contacting lower 4 incisors and extends labially and lingually to patients range of motion
- * No posterior tooth contact.

We take great pride in not only meeting, but exceeding the expectations of our clients. We hope you will agree when you try our latest dental orthotics. If you have any questions about orthotics, please contact us.