



**Klausz Dental
Laboratories Ltd.**

Working **harder and smarter** for your practice!

www.klauszdentallab.com • klauszdl@uniserve.com • Tel: 416-630-6262 • 1-800-410-4392
40 Kodiak Crescent • Units 11-12 • Toronto, ON • Canada M3J 3G5

An O.I.C.S. Specialty Laboratory
Occlusion - Implants - Cosmetic - Sleep

A Simple, Strong, Comfortable, Adjustable, & Low Cost Sleep Apnea Device

Snoring can be a real problem, not only for the person who snores, but also for their partner and other family members. However, there is a dentist-prescribed solution for a more restful night's sleep: QuietNite. This affordable device is flexible, thin and comfortable for your patients, and it exhibits clinical success in mitigating or even preventing the disruptive, unhealthy effects of snoring and sleep apnea.

QuietNite is a mandibular re-positioning and advancement device which is intended for the treatment of snoring and mild to moderate obstructive sleep apnea in persons 18 years of age or older.

QuietNite consists of two separate thermal forming trays, which fit over the upper and lower arches. These trays are comprised of a soft inner layer with a hard outer layer that is durable and BPA-free. The trays are connected by two lateral rotating connectors, which gently pull the lower jaw forward. The connectors are easily interchangeable with longer or shorter connectors. They also exhibit extraordinary strength when subjected to lateral forces such as bruxism, significantly reducing breakage/failures and improving reliability. Patients can adjust the connectors to meet their comfort level at home.

QuietNite devices do not interfere with breathing through the mouth. Small movements of the jaw (TMJ) are possible while wearing a QuietNite device. This movement potentially helps minimize stiffness in the jaw joints in the morning.

Research has shown that custom fabricated dental devices worn at night that move the lower jaw into a forward position increase the three dimensional space in the airway tube, which reduces air velocity and soft tissue vibration. By increasing the volumetric capacity of the airway and preventing soft tissue vibrations, snoring can be reduced or even eliminated. QuietNite devices are a simple, cost effective way to help treat your patients.

Quiet
NITE



Every day we strive to exceed your expectations and being able to provide you with another great product is just another way we are, "Working Harder and Smarter for your Practice!"