



**Klausz Dental
Laboratories Ltd.**

Working **harder and smarter** for your practice!

www.klauszdentallab.com • klauszdl@uniserve.com • Tel: 416-630-6262 • 1-800-410-4392
40 Kodiak Crescent • Units 11-12 • Toronto, ON • Canada M3J 3G5

An O.I.C.S. Specialty Laboratory
Occlusion - Implants - Cosmetics - Sleep

Splint Therapy Just Got A Lot More Fun!!!

Our Thermoplastics are now available in a variety of colours, designed to be as individual as your patients are.



We are proud to be able to offer our clients a variety of options when it comes to splint therapy, and now all our Thermoplastic Splints can be tinted in 8 additional colours to complement our original crystal clear shade.

Clear Pink Blue Orange Yellow Violet Green Black Red

The Thermoplastic Splint is made from a unique thermo-plastic material and was developed for greater patient comfort. The advanced properties of this material allow it to flex when softened under warm water and thus conforms to the contours of the patient's teeth when inserted.

The Hard Occlusal/Thermoplastic Interior Splint is by far the most popular splint and was developed for those heavy bruxers and clenchers. Its unique combination of a tough impact resistant occlusal surface with its pliable tooth and tissue side offers your patients the best of both worlds: lasting durability and custom comfort.

The Thermoplastic Anterior Deprogrammer is designed to prevent posterior tooth contact. It is more than a standard bite splint. It protects the teeth by calming the temporalis muscles, and in many cases can reduce or eliminate muscle pain and inflammation associated with bruxing and clenching.

Whatever your requirements, we have the splint for you! Now In Colour!