



**Klausz Dental  
Laboratories Ltd.**

Working **harder and smarter** for your practice!

www.klauszdentallab.com • klauszdl@uniserve.com • Tel: 416-630-6262 • 1-800-410-4392  
40 Kodiak Crescent • Units 11-12 • Toronto, ON • Canada M3J 3G5

## **An O.I.C.S. Specialty Laboratory** *Occlusion - Implants - Cosmetic - Sleep*

# MiniMarvel Appliance Developed to Reduce Daytime Clenching and Grinding

The MiniMarvel Appliance is a discreet custom fabricated splint for your patients who are looking for ways to stop daytime clenching and grinding.



Experiences involving physical exertion, stress, anger, and intense concentration all have a higher potential for clenching damage.

The MiniMarvel Appliance design is light, thin, and very comfortable. This appliance is so discreet; it doesn't even affect speech. The appliance easily trains the wearer to keep his or her teeth apart.

### **Exceptional Benefits of the MiniMarvel Appliance**



- Recommended for Daytime Wear – The appliance is almost invisible when worn.
- No Adjustments – It seats with ease, with no time wasted on grinding in the occlusion.
- Improved Compliance and Performance – It is extremely comfortable and very stable, due to its use of tissue coverage to create a Hydroseal.

The secret to its comfort lies in the strategic placement of two discluding elements which act as rigid therapeutic bands. These are positioned to occlude with the maxillary canines making sure that the posterior teeth are separated. These elements deprogram the muscles of mastication by minimizing bite force, thus preventing the patient from completely occluding, and helping maintain a comfortable position. Since there are just two therapeutic bands as points of occlusal contact, there's no need for any adjustments.



Different from other splints, MiniMarvel Appliances will remind the patient when they are biting too hard, in an effort to change the behavior.

Simple and effective, the MiniMarvel Appliance may just be what you and your patients are looking for. Looking for innovative ideas, is just another way we are, "Working Harder and Smarter for your Practice!"