



**Klausz Dental
Laboratories Ltd.**

Working **harder and smarter** for your practice!

www.klauszdentallab.com • klauszdl@uniserve.com • Tel: 416-630-6262 • 1-800-410-4392
40 Kodiak Crescent • Units 11-12 • Toronto, ON • Canada M3J 3G5

FALL 2023

INSIDE THIS ISSUE

- 1** We've Grown into an O.I.C.S. Specialty Laboratory
- 1** Are You Looking For More?
- 2** Oral Appliance Care
- 2** Splint Therapy Just Got a Lot More Fun!
- 3** Let's Stop Wasting Chair Time With Our Exclusive Bite Equilibration Technique
- 4** O2Vent Optima – Airway Management System
- 4** Aqualizer 7th Annual November Sale

"We Build Advanced Occlusal Concepts Into Everything We Do!"

We've grown into an O.I.C.S. Specialty Laboratory!

At **Klausz Dental Laboratories**, we are committed to providing our clients with a full range of services, restorations, and appliances to meet all aspects of clinical dentistry. To meet your expectations we have grown into an **Occlusal-Implant-Cosmetic-Sleep (O.I.C.S.) Specialty Laboratory** and we have solidified our position as a leader in these fields through various certifications.

With over 55 years of experience in the dental technology field, we build advanced occlusal concepts into every restoration and appliance we custom fabricate. It's with a greater understanding of how teeth, muscles and joints work together that we can create restorations that not only look lifelike and aesthetic, but can last longer and function better. Allow us to apply the knowledge we've accumulated into your next restoration or appliance. See why we've become an **O.I.C.S. Specialty Laboratory**.

Are You Looking For More?

While most labs are busy trying to cut costs and find shortcuts, we're busy looking for ways to improve on our quality and service.

With our **Exclusive Bite Equilibration Technique**, applied to all of our fixed restorations, we have been able to successfully eliminate approximately 95% of all chair side adjustments.

Our own **Modified Denture Technique** reduces the need for several appointments and reduces the fees for complete denture patients.

And now, with all our removable appliances, we include a free sample of our private line of **Oral Appliance Care Dental Appliance Cleaner**, to ensure the appliance we fabricate today can be properly cared for.

We're looking to build lasting relationships with our clients. When we say, "Working Harder and Smarter For Your Practice", we mean that we're finding solutions for each of our client's specific needs and challenges, whatever they may be.



Ronald Klausz, R.D.T.

Are you looking for more? If so, why not give me a call!

Helping To Care For The Products We Create With Innovative Cleaning And Care Solutions

New Name - Same Great Product



We've always been dedicated to providing our clients, and their patients, with superior dental appliances. We are equally committed to ensuring that these appliances can be properly cared for. For years, we've been including a sample of our exclusively formulated "Spray'n Fresh Dental Appliance Cleaner" with all our splints, night guards, sports guards, orthotics, sleep appliances, and retainers. We've rebranded this product to improve its appeal with a fresh new look and improved mint scent. It's the same great product, just with a new look.

Oral Appliance Care Dental Appliance Cleaner can be easily reordered from our newly created website at:

www.oralappliancecare.com

Splint Therapy Just Got a Lot More Fun!

Our Thermoplastics are now available in a variety of colours, designed to be as individual as your patients are.

We are proud to be able to offer our clients a variety of options when it comes to splint therapy, and now all our Thermoplastic Splints can be tinted in 8 additional colours to complement our original crystal clear shade.



Clear Pink Blue Orange Yellow Violet Green Black Red

The **Thermoplastic Splint** is made from a unique thermo-plastic material and was developed for greater patient comfort. The advanced properties of this material allow it to flex when softened under warm water and thus conforms to the contours of the patient's teeth when inserted.

The **Hard Occlusal/Thermoplastic Interior Splint** is by far the most popular splint and was developed for those heavy bruxers and clenchers. Its unique combination of a tough impact resistant occlusal surface with its pliable tooth and tissue side offers your patients the best of both worlds: lasting durability and custom comfort.

Whatever your requirements, we have the splint for you! Now In Colour!

Your portal to all that's new and exciting in dental technology:



www.klauszdentallab.com



Let's Stop Wasting Chair Time With Our Exclusive Bite Equilibration Technique!

When was the last time you inserted a crown with very little occlusal adjustments? If your answer is "Never", then you need to experience the best of crown and bridge, combined with our **Exclusive Bite Equilibration Technique**, which was developed to eliminate most chair side occlusal adjustments.

In a poll of dental offices, the number one complaint with fixed restorations, be it porcelain fused to metal, all ceramics, or polymer glass is the excessive time loss associated with adjusting the occlusion or even having to grind in the bites.

We've found a solution that has proven itself over time. With our **Exclusive Bite Equilibration Technique**, we have been able to successfully eliminate approximately 95% of all chair side occlusal adjustments.

The technique consists of eliminating all bubbles from the occlusal surfaces and with fine articulating paper, looking for high spots in occlusion. We then eliminate these high spots until more stable contact markings appear. Now the mounted models more accurately represent the patient's occlusion.




If you are interested in finding out how you can benefit from this technique, just give us a call and we will be pleased to discuss it further. Better yet, send us a case and we will show it to you. The sooner you call, the sooner you can stop adjusting bites and start working more efficiently.

Would You Like to Reduce Your Chair
Side Occlusal Adjustments by 95%?

If so, be sure to read Ron Klausz's article
which was published in the
Oral Health Journal

<https://www.klauszdentallab.com/published-articles/>





New from Klausz
Dental
Laboratories

O2Vent Optima, Airway Management System

O2Vent Optima

The O2Vent Optima appliances are the latest offering in the O2Vent oral appliance range. O2Vent Optima is 3D printed with durable nylon. The lower jaw can be moved forward to a comfortable range while also allowing the integrated airway channel to bypass nasal congestion and other areas of airway collapse, such as from the tongue, soft palate and lateral walls.



O2Vent Optima Mini

The O2Vent Optima Mini is a discreet and low-profile oral device that manages the upper airway for the treatment of mild to moderate obstructive sleep apnea (OSA) and snoring by incorporating the proprietary Oventus Airway Technology. The O2Vent Optima Mini's design (without the airway extension) may benefit patients requiring enhanced lip-seal, to avoid excess saliva build-up and/or dry lips.



Oventus ExVent

ExVent is an optional accessory which is inserted into the front oral airway of the O2Vent Optima appliance. It creates resistance in the upper airway, providing additional airway support whilst prolonging exhalation.

How ExVent works

ExVent - Breathing in

ExVent - Breathing out



*ExVent is not a stand-alone device and cannot be used with other mandibular advancement splints (MAS).

Always looking for new and innovative appliances to help your patients sleep better is just another way we are, "Working Harder and Smarter for your Practice"

It's back again - Our 7th Annual Aqualizer Sale is once again being held in November!



Laboratory Sponsor

www.SleepDisordersDentistry.com

Working with others to provide our clients with the very best in Continuing Education is just another way we are, "Working Harder and Smarter for your Practice!"



Laboratory Sponsor

www.cacsdd.org